Vertrieb/Distributeur:
Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo - Germany
Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-289
info@kometdental.de
www.kometdental.de

Hersteller/Manufacturer:
Elaboro GmbH
Hagenower Straße 73
19061 Schwerin - Germany
Telefon +49 (0) 385 3993-300
Telefax +49 (0) 385 3993-302
info@elaboro.de
www.elaboro.de

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Second firing (to achieve customized results)

| Closure time | Stand-by temperature | Heating rate | Firing temperature | Hold time | Heating rate | Firing temperature | Hold time | Slow cooling | Cooling rate | Vacuum 1 |
|--------------|----------------------|--------------|--------------------|-----------|--------------|--------------------|-----------|--------------|--|----------|
| min. | °C | °C/min. | °C/min. | min. | °C/min. | °C | min. | °C/min. | According to the size of the framework | 400°C |
| 1 | 400 | - | - | - | 50 | 920°C | 2 | Ja | 20 ... 80°C/min. | max. 30% |

Alternative firing

| Closure time | Stand-by temperature | Heating rate | Firing temperature | Hold time | Slow cooling | Cooling rate | Vacuum 1 |
|--------------|----------------------|--------------|--------------------|-----------|--------------|--|----------|
| min. | °C | °C/min. | °C | min. | °C/min. | According to the size of the framework | off |
| 1 | 400 | 50 | 920°C | 5 | Ja | 20 ... 80°C/min. | off |

CeraFusion evo

Typing, classification and identification as per DIN EN ISO 6872; Dental ceramic type I, class 1a, classification according to Directive 93/42/EEC, annex IX, rule 6: Medical device IIa

1. Scope of application

CeraFusion evo is a ready-to-use and easy-to-apply lithium silicate based glass ceramic for monolithic restorations made of zirconium oxide. To optimally coat the prosthesis with CeraFusion evo, no special surface preparation, adhesion promoter, solvents or preparatory firing are required. The parts made of zirconium oxide should be clean and free of dust and grease to ensure an even result. CeraFusion evo is suitable for all dental zirconium oxides and many high-fusion framework and veneer ceramics. Best results are achieved by finishing highly translucent and colored ceramic frameworks made of zirconium oxide or pressed ceramics.

2. Working environment

CeraFusion evo spray bottles should be stored at room temperature. CeraFusion evo should be used at moderate ambient temperatures of 15-35°C. Very hot or low temperatures should be avoided because they have a negative effect on the result. Only spray in well ventilated rooms, use suitable exhaust systems and wear a dust mask to prevent inhalation of the fine spray mist. Good lighting is essential for checking that the prosthesis is evenly coated without gaps.

3. General instructions for handling

- CeraFusion evo is intended solely for use by trained and qualified personal in dental laboratories.
- Only use in well ventilated rooms
- Do not breathe in the spray mist
- Always wear a dust mask and use a suitable exhaust system
- Do not use inside the mouth
- Follow the instructions on the material safety data sheet.
- The aerosol can is under pressure and must be protected from solar radiation and temperatures above 50°C
- Keep away from ignition sources.
- Do not smoke.
- Do not spray into open flames (Bunsen burner, for example) or any incandescent material.
- Keep the spray out of the reach of children.
- Do not open aerosol can with force or incinerate after use.
- Always completely empty the aerosol can.

4. Shaking of the aerosol can

Step 1: To activate the components of the spray, shake can vigorously before first applying the spray head to ensure that the ceramic particles mix with the mixing liquid in the spray can. Although you can hear the mixing balls after a few seconds, continue shaking the can vigorously for 1 minute.

Step 2: Place the spray head onto the spray lance and briefly shake again. The spray can is now ready for use. After every break, shake the can again performing circular movements. This procedure is absolutely essential and has to be followed without fail for optimum results and to prevent failure of the nozzles, the riser tube and the valve system.

5. Preparation of the zirconium oxide restoration

To achieve an optimum surface coating, the ready sintered and fitted restoration made of zirconium oxide has to be dry, clean and free of dust and grease. Only a small amount of CeraFusion evo needs to be applied to coat surfaces. Functional and aesthetic characteristics such as the occlusal surface, contact points, color gradients and effects should already be considered prior to sintering the zirconium oxide, especially in case of fully contoured monolithic restorations. The occlusal relief, the crown margins and the margin fitting are hardly affected by one application of CeraFusion evo. We therefore recommend not to apply CeraFusion evo until after the denture has been fitted and any corrections have been made.

6. Further indications

Surface treatment of veneering framework:
A thin layer of CeraFusion evo is applied onto the surface of the veneers and fired according to the instructions for diffusion firing.

Surface treatment of adhesive constructions:

A thin layer of CeraFusion evo is applied on the inner surface of the crown or on the wing area of a Maryland bridge and fire according to the instructions for diffusion firing. Coat and fire the inner surface first, and then the outer surface. Like this, your zirconium restoration is optimally prepared for adhesive fixation.

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© Gebr. Brasseler GmbH & Co. KG