Vertrieb/Importeur:
Komet Dental GmbH & Co. KG
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo - Germany
Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-289
info@kometdental.de
www.kometdental.de

Hersteller/Hersteller:
Elaboro GmbH
Hagenower Straße 73
19061 Schwerin - Germany
Telefon +49 (0) 385 3993-300
Telefax +49 (0) 385 3993-302
info@elaboro.de
www.elaboro.de

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Second firing (to achieve customized results)

Closure time	Stand-by temperature	Heating rate	Firing temperature	Hold time	Heating rate	Firing temperature	Hold time	Slow cooling	Cooling rate	Vacuum 1
min.	°C	°C/min.	°C/min.	min.	°C/min.	°C	min.	°C/min.	According to the size of the framework	400°C
1	400	-	-	-	50	920°C	2	Ja	20 ... 80°C/min.	max. 30%

Alternative firing

Closure time	Stand-by temperature	Heating rate	Firing temperature	Hold time	Slow cooling	Cooling rate	Vacuum 1
min.	°C	°C/min.	°C	min.	°C/min.	According to the size of the framework	off
1	400	50	920°C	5	Ja	20 ... 80°C/min.	off

CeraFusion evo

Typing, classification and identification as per DIN EN ISO 6872: Dental ceramic type I, class 1a, classification according to Directive 93/42/EEC, annex IX, rule 6: Medical device IIa

1. Scope of application

CeraFusion evo is a ready-to-use and easy-to-apply lithium silicate based glass ceramic for monolithic restorations made of zirconium oxide. To optimally coat the prosthesis with CeraFusion evo, no special surface preparation, adhesion promoter, solvents or preparatory firing are required. The parts made of zirconium oxide should be clean and free of dust and grease to ensure an even result. CeraFusion evo is suitable for all dental zirconium oxides and many high-fusion framework and veneer ceramics. Best results are achieved by finishing highly translucent and colored ceramic frameworks made of zirconium oxide or pressed ceramics.

2. Working environment

CeraFusion evo spray bottles should be stored at room temperature. CeraFusion evo should be used at moderate ambient temperatures of 15-35°C. Very hot or low temperatures should be avoided because they have a negative effect on the result. Only spray in well ventilated rooms, use suitable exhaust systems and wear a dust mask to prevent inhalation of the fine spray mist. Good lighting is essential for checking that the prosthesis is evenly coated without gaps.

3. General instructions for handling

CeraFusion evo is intended solely for use by trained and qualified personal in dental laboratories.
- Only use in well ventilated rooms
- Do not breathe in the spray mist
- Always wear a dust mask and use a suitable exhaust system
- Do not use inside the mouth
- Follow the instructions on the material safety data sheet.
- The aerosol can is under pressure and must be protected from solar radiation and temperatures above 50°C
- Keep away from ignition sources.
- Do not smoke.
- Do not spray into open flames (Bunsen burner, for example) or any incandescent material.
- Keep the spray out of the reach of children.
- Do not open aerosol can with force or incinerate after use.
- Always completely empty the aerosol can.

4. Shaking of the aerosol can

Step 1: To activate the components of the spray, shake can vigorously before first applying the spray head to ensure that the ceramic particles mix with the mixing liquid in a spray can. Although you can hear the mixing balls after a few seconds, continue shaking the can vigorously for 1 minute.

Step 2: Place the spray head onto the spray lance and briefly shake again. The spray can is now ready for use. After every break, shake the can again performing circular movements. This procedure is absolutely essential and has to be followed without fail for optimum results and to prevent failure of the nozzles, the riser tube and the valve system.

5. Preparation of the zirconium oxide restoration

To achieve an optimum surface coating, the ready sintered and fitted restoration made of zirconium oxide has to be dry, clean and free of dust and grease. Only a small amount of CeraFusion evo needs to be applied to coat surfaces. Functional and aesthetic characteristics such as the occlusal surface, contact points, color gradients and effects should already be considered prior to sintering the zirconium oxide, especially in case of fully contoured monolithic restorations. The occlusal relief, the crown margins and the margin fitting are hardly affected by one application of CeraFusion evo. We therefore recommend not to apply CeraFusion evo until after the denture has been fitted and any corrections have been made.

6. Further indications

Surface treatment of veneering framework:
A thin layer of CeraFusion evo is applied onto the surface of the veneers and fired according to the instructions for diffusion firing.

Surface treatment of adhesive constructions:
A thin layer of CeraFusion evo is applied on the inner surface of the crown or on the wing area of a Maryland bridge and fire according to the instructions for diffusion firing. Coat and fire the inner surface first, and then the outer surface. Like this, your zirconium restoration is optimally prepared for adhesive fixation.

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Stains and veneers:
CeraFusion evo is a universal, colorless veneer that is applied in a thin layer. It can be used on its own or in combination with high-fusion stains and ceramic masses. If the restoration is to be customized with stains, this should be applied directly after the first firing (diffusion firing). Once the stain has dried, another thin layer of CeraFusion evo is applied and the restoration is fired for the second time. After the second firing, the surface should be shiny, homogenous and free of pores.

7. Spraying

Prior to first use, spray a small amount onto a piece of glass to test the effect. Make sure that the correct distance to the object, i.e. approx. 10-20 cm is observed. For best results, spray on in small, short bursts. Wet patches or knobs occur if the product has been sprayed on from too close a distance. Uneven deposits of powder are signs of excessive spraying or insufficiently mixed powder.

Layer thickness:

Do not spray excessively—always make sure that the zirconium surface still shines through the white powder. The carrier fluid evaporates completely after a few seconds and leaves a thin layer of powder which adheres well to the surface and does not run.

Spraying techniques:

Hold the spray bottle as vertically as possible. Unlike during "varnishing", it is advisable to apply CeraFusion evo in small short spray bursts. Like this, only small amounts of powder are dispersed and the nozzle system cleans itself at the same time.

8. Cleaning

Clean the spray head immediately after use (for example, with a steam jet air ejector or ultrasonic bath) and dry with compressed air.

9. Ceramic firing

Carry out the ceramic firing according to the instructions (firing parameters) for CeraFusion evo.

Diffusion firing:

If the product is applied correctly, a resilient, smooth and homogenous surface of outstanding quality and a perfect bond with the zirconium oxide is achieved after just one firing.

Firing to achieve customized results:

During the second firing, the hold times can be shortened to reduce the thermal loads exerted on the stains.

Alternative firing recommendation:

We have added a simplified firing program with a certain heating rate for older ceramic furnaces. This is appreciated by many customers.

We suggest that you define the parameters that a best suited to your own personal requirements with an old crown that is no longer used.

10. Troubleshooting

Incorrectly or inadvertently sprayed parts can easily be cleaned with water or steam. Dry the steamed surface with compressed air and spray on another layer immediately afterwards. Small pores in the surface after firing often occur when the layer of powder applied is too thin. In this case, simply spray another layer of CeraFusion evo onto the prosthesis without any previous surface treatment and repeat firing. Check the firing parameters, if necessary.

Please observe the instructions on the label of the spray can:

Flammable aerosol. Pressurised container; may burst if heated. Do not spray on an open flame or other ignition source. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Do not pierce or burn, even after use. Protect from sunlight. Do not expose to temperatures exceeding 50°C. Keep out of reach of children. May irritate mucosa of nose and throat.

The firing parameters specified in the instructions for use have to be observed without fail.

11. Instructions for the user

Surfaces treated with CeraFusion evo must not be blasted with the powdery blasting agent during professional scaling and polishing. The roughening of the glass ceramic is irreversible and can increase the susceptibility of plaque deposits settling on the surfaces.

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© Gebr. Brasseler GmbH & Co. KG