

Rotadent

ROTATING INSTRUMENTS



Endodontické nástroje - VDW

Pouze pro potřeby firmy Rotadent spol.s r.o.

We bring good *shapes* to life



RECIPROC® blue



A great file. Even better.

- Same cutting efficiency
- Increased flexibility and safety
- Same simplicity of use


Simple, but effective: The new RECIPROC® blue file generation combines the ease of the original RECIPROC® one file endo concept with enhanced safety in root canal preparation for patients. An innovative heat treatment makes RECIPROC® blue particularly flexible to ensure a smoother and safer progression in the canal and gives it in addition its characteristic blue color.



RECIPROC® blue Single Sizes

	Blister pack of 6 instruments	21 mm	25 mm	31 mm
●	R25	(6632) V04 0252 021 025	(6636) V04 0252 025 025	(6640) V04 0252 031 025
●	R40	(6633) V04 0252 021 040	(6637) V04 0252 025 040	(6641) V04 0252 031 040
●	R50	(6634) V04 0252 021 050	(6638) V04 0252 025 050	(6642) V04 0252 031 050

STERILE



RECIPROC® blue Assortments

	Blister pack of 6 instruments	21 mm	25 mm
	3 × R40, 3 × R50	(6635) V04 0052 021 200	(6639) V04 0052 025 200

STERILE



RECIPROC®

System designed specifically for root canal preparation in reciprocation

- Root canal preparation with only one instrument
- Retreatment of gutta-percha filling material and carrier based obturators with R25
- Single-use convenience; sterile instruments, individually packed
- Three instrument sizes for narrow, medium and wide canals

RECIPROC® Single Sizes

STERILE

	Blister pack of 6 instruments	21 mm	25 mm	31 mm
●	R25	(6620) V04 0212 021 025	(6624) V04 0212 025 025	(6628) V04 0212 031 025
●	R40	(6621) V04 0212 021 040	(6625) V04 0212 025 040	(6629) V04 0212 031 040
●	R50	(6622) V04 0212 021 050	(6626) V04 0212 025 050	(6630) V04 0212 031 050



RECIPROC® Assortments

STERILE

	Blister pack of 6 instruments	21 mm	25 mm	31 mm
	3 × R40, 3 × R50	(6623) V04 0212 021 233	(6627) V04 0212 025 233	(6631) V04 0212 031 233



R-PILOT™



První reciproční nástroj pro glide path od firmy VDW


Safer and faster glide path management** in reciprocating motion

- Better respect of the canal anatomy** for more confidence in your treatment

With R-PILOT™, there is finally a safe and easy-to-use solution for cases where a glide path is needed. R-PILOT™ is designed to be used in reciprocating motion with a designated drive system with original VDW RECIPROC® settings.

STERILE

	Blister pack of 6 instruments	21 mm	25 mm	31 mm
●	12.5/.04	(8043) V04 0217 021 012	(8044) V04 0217 025 012	(8845) V04 0217 031 012



M_{two}® Retreatment Files

STERILE

For efficient removal of filling material

Suitable for narrow, medium and large root canals

- Cutting tip
- Working part 16 mm

	Size/Taper	21 mm
○	R 15/.05	(6696) V04 0232 021 015
●	R 25/.05	(6697) V04 0232 021 025
	3 × R 15/.05, R 25/.05 each	(6698) V04 0232 021 202



VDW®STERILE C-PILOT® Files

STERILE

STERILE

For severely curved or calcified root canals

- Stainless steel
- High bending strength due to a special thermal hardening process
- Enables creation of a smooth glide path prior to every rotary preparation
- The inactive pilot tip conducts the instrument safely along the canal
- Blister pack of 6 pieces

	Size	19 mm	21 mm	25 mm
○	06	(6788) V04 0368 019 006	(6794) V04 0368 021 006	(6800) V04 0368 025 006
○	08	(6789) V04 0368 019 008	(6795) V04 0368 021 008	(6801) V04 0368 025 008
○	10	(6790) V04 0368 019 010	(6796) V04 0368 021 010	(6802) V04 0368 025 010
	06-10	(6793) V04 0368 019 210	(6799) V04 0368 021 210	(6805) V04 0368 025 210
○	12.5	(6791) V04 0368 019 012	(6797) V04 0368 021 012	(6803) V04 0368 025 012
○	15	(6792) V04 0368 019 015	(6798) V04 0368 021 015	(6804) V04 0368 025 015



FlexMaster®

The established crown-down NiTi system

- Three sequences to choose from for every type of canal
- Safe use due to optimal stability of file even with low tapered files (2 %)
- Special K-type cutting edges for efficient removal of debris and dentine chips

FlexMaster® kónus .02 blister, baleno po 6 kusech

	Size/Taper	21 mm	25 mm
○	15/.02	(8156) V04 0341 021 015	(8163) V04 0341 025 015
●	20/.02	(8157) V04 0341 021 020	(8164) V04 0341 025 020
●	25/.02	(8158) V04 0341 021 025	(8165) V04 0341 025 025
●	30/.02	(8159) V04 0341 021 030	(8166) V04 0341 025 030
●	35/.02	(8160) V04 0341 021 035	(8167) V04 0341 025 035
●	40/.02	-	(8168) V04 0341 025 040
○	45/.02	-	(8169) V04 0341 025 045
●	50/.02	-	(8170) V04 0341 025 050
●	60/.02	-	(8171) V04 0341 025 060




STERILE



Mtwo® is licensed by SWEDEN & MARTINA

The efficient NiTi system

- Simple: One sequence for all canals
- Single Length: All instruments are brought to full working length
- Simultaneous Shaping: Mtwo® cuts automatically while advancing apically and also laterally when using a brushing file movement



Basic Sequence Blister pack-baleno po 6 nástrojích



		Mtwo® pracovní délka 16 mm		Mtwo® pracovní délka 21 mm	
Size/Taper		21 mm	25 mm	25 mm	31 mm
10/.04	(6644) V04 0234 021 010 (6664) V04 0234 025 010			(6678) V04 1234 025 010	(6682) V04 1234 031 010
15/.05	(6645) V04 0235 021 015 (6665) V04 0235 025 015			(6679) V04 1235 025 015	(6683) V04 1235 031 015
20/.06	(6646) V04 0236 021 020 (6666) V04 0236 025 020			(6680) V04 1236 025 020	(6684) V04 1236 031 020
25/.06	(6652) V04 0236 021 025 (6667) V04 0236 025 025			(6681) V04 1236 025 025	(6685) V04 1236 031 025



Opracování větších kanálků

Blister pack-baleno po 6 nástrojích



		Mtwo® pracovní délka 16 mm		Mtwo® pracovní délka 21 mm	
Size/Taper		21 mm	25 mm	31 mm	
30/.05	(6653) V04 0235 021 030 (6668) V04 0235 025 030			(6686)	V04 1235 031 030
35/.04	(6655) V04 0234 021 035 (6669) V04 0234 025 035			(6687)	V04 1234 031 035
40/.04	(6656) V04 0234 021 040 (6670) V04 0234 025 040			(6688)	V04 1234 031 040
45/.04	(6657) V04 0234 021 045 (6671) V04 0234 025 045			(6689)	V04 1234 031 045
50/.04	(6658) V04 0234 021 050 (6672) V04 0234 025 050			(6690)	V04 1234 031 050
60/.04	(6659) V04 0234 021 060 (6673) V04 0234 025 060			(6691)	V04 1234 031 060
		Mtwo® pracovní délka 16 mm		Mtwo® pracovní délka 21 mm	
Size/Taper		21 mm	25 mm	31 mm	
25/.07	(6660) V04 0237 021 025 (6674) V04 0237 025 025			(6692)	V04 1237 031 025
30/.06	(6661) V04 0236 021 030 (6675) V04 0236 025 030			(6693)	V04 1236 031 030
35/.06	(6662) V04 0236 021 035 (6676) V04 0236 025 035			(6694)	V04 1236 031 035
40/.06	(6663) V04 0236 021 040 (6677) V04 0236 025 040			(66995)	V04 1236 031 040



FlexMaster®

The established crown-down NiTi system

- Three sequences to choose from for every type of canal
- Safe use due to optimal stability of file even with low tapered files (2 %)
- Special K-type cutting edges for efficient removal of debris and dentine chips

FlexMaster® kónus .04 blister, baleno po 6 kusech

	Size/Taper	21 mm	25 mm	31 mm
○	15/.04	-	(8179) V04 0342 025 015	-
●	20/.04	-	(8180) V04 0342 025 020	(8185) V04 0342 031 020
●	25/.04	(8175) V04 0342 021 025	(8181) V04 0342 025 025	(8186) V04 0342 031 025
●	30/.04	-	(8182) V04 0342 025 030	(8187) V04 0342 031 030
●	35/.04	-	(8183) V04 0342 025 035	-
●	40/.04	-	(8184) V04 0342 025 040	-



FlexMaster® kónus .06 blister, baleno po 6 kusech

	Size/Taper	25 mm
○	15/.06	
●	20/.06	(8197) V04 0343 025 020
●	25/.06	(8198) V04 0343 025 025
●	30/.06	(8199) V04 0343 025 030
●	35/.06	(8200) V04 0343 025 035
●	40/.06	(8201) V04 0343 025 040



FlexMaster® IntroFile

- 11 % taper
- 19 mm pracovní délka, 9 mm aktivní pracovní délka
- ISO size 22 at the tip
- Blister, baleno po 6 kusech

(8219) V04 0357 019 000




EDDY®

Innovative sonic powered irrigation – safe, effective and proven

- Sonic powered polyamide tip
- Powerful activation at max. 6,000 Hz via Airscaler
- Reliable removal of residual tissue and dentine chips
- Three-dimensional movement for irrigant activation
- Many clinical studies substantiate the effectiveness of EDDY®

Cavitation

Due to the fast movement of the irrigating tip in the irrigant small gas bubbles form and then implode. The energy released intensifies the cleaning effect of the irrigant.



Acoustic Streaming

High frequency acoustic oscillation creates swirling eddies in the irrigant which strengthens the cleaning effect on the canal walls.

› EDDY® _____ (8064) V04 1441 000 000

› EDDY® Starter Kit _____ (8065) V04 1443 000 000

- 1 W&H Airscaler Handpiece, Type Proxeo ZA-55-LM (Multiflex®) or ZA-55-LS (Sirona®) with light and power setting
- 4 balení po 10 EDDY® špičkách, sterilní

(Multiflex®)

(Sirona®)

› EDDY® 100 Plus Kit _____ (8067) V04 6034 000 000

- 60 blisterů baleno po 2 EDDY® špičkách, sterilní



RECIPROC® blue Gutta-Percha

Suitable for the RECIPROC® blue NiTi system

- Multi-variable taper for increased precision and consistency to match corresponding file system
- Latex-free
- Heat Extension < 4 mm
- Flip pack of 60 pieces

	Size	28 mm
●	R25	(8068) 709 V04 0258 028 025
●	R40	(8069) 709 V04 0258 028 040
●	R50	(8070) 709 V04 0258 028 050
36 × R25, 12 × R40, 12 × R50		(8071) 709 V04 0258 028 237



RECIPROC® Gutta-Percha

Suitable for the RECIPROC® NiTi system Branded system solution for simple integration

- Slide cassette of 60 pieces

	Size	28 mm
●	R25	(8072) 710 V04 0214 028 025
●	R40	(8073) 710 V04 0214 028 040
●	R50	(8074) 710 V04 0214 028 050
40 × R25, 10 × R40, 10 × R50		(8075) 710 V04 0214 028 237



Mtwo® Gutta-Percha

Suitable for the Mtwo® NiTi system

- α-phase with low melting point for cold and warm obturation techniques
- Slide cassette of 60 pieces

	Size/Taper	28 mm
●	35/.04	(8076) 711 V04 0220 428 035
●	40/.04	(8077) 711 V04 0220 428 040
●	45/.04	(8078) 711 V04 0220 428 045
●	50/.04	(8079) 711 V04 0220 428 050
●	60/.04	(8080) 711 V04 0220 428 060
●	30/.05	(8081) 711 V04 0220 528 030
●	25/.06	(8082) 711 V04 0220 628 025
●	30/.06	(8083) 711 V04 0220 628 030
●	35/.06	(8084) 711 V04 0220 628 035
●	40/.06	(8085) 711 V04 0220 628 040
●	25/.07	(8086) 711 V04 0220 728 025
20 × 25/.06, 10 × 30/.05, 35/.04, 40/.04, 25/.07 each		(8087) 711 V04 0220 028 457



BeeFill® 2in1

BeeFill® Heat Pluggers

- For Downpack and searing off gutta-percha
- Continuous Wave or Schilder technique
- Fast heating and cooling down
- Autoclavable

Heat plugger standard 40/.03 (8202) V04 0547 000 040

Heat plugger soft 50/.05 (8203) V04 0546 000 050

Heat plugger soft 60/.06 (8204) V04 0546 000 060



BeeFill® Cartridges

- Injection of warm gutta-percha
- Flexible and extra long canula – also for curved root canals
- One cartridge fills 4-6 root canals
- Matching cartridge for each canal anatomy:
25G = Ø 0.45 mm | 23G = Ø 0.60 mm | 20G = Ø 0.80 mm

Size

20G / 0.80 mm (8205) V04 0532 053 080

23G / 0.60 mm (8206) V04 0532 053 060

25G / 0.45 mm (8207) V04 0532 053 045



Příslušenství

› Minifix Measuring Gauge

For precise adjustment of the working length of endo instruments, (8046) 722 V04 0179 gutta-percha and paper points



› Čistící stojánek

For quick chairside storing and cleaning of root canal instruments during treatment

- Foam disc autoclavable for one single use
- Čistící stojánek je sterilizovatelný, ale bez molitanových vložek do teploty 200°C (8047) 723 V04 0495
- Bez nástrojů



› Náplně do čistícího stojánku

- Foam discs autoclavable for one single use (8048) 723 V04 0496
- 55 kusů



RECIPROC® *blue* Paper Points

STERILE

Suitable for the RECIPROC® *blue* NiTi system

- Fully saturated points retain good material consistency
- Sterilized points in blister cells to support infection control
- Part of a system based solution for RECIPROC® *blue*
- 180 pieces

	Size/Taper	29 mm
●	R25	(8088) V04 0259 029 025
●	R40	(8089) V04 0259 029 040
●	R50	(8090) V04 0259 029 050
	60 × R25, 60 × R40, 60 × R50	(8091) V04 0259 029 237



RECIPROC® Paper Points

STERILE

Suitable for the RECIPROC® NiTi system

- Fully saturated points retain good material consistency
- Sterilised points in blister cells to support infection control
- 144 pieces: 4 papírové čepy v jednom blistru

	Size/Taper	29 mm
●	R25	(8092) V04 0216 029 025
●	R40	(8093) V04 0216 029 040
●	R50	(8094) V04 0216 029 050
	96 × R25, 24 × R40, 24 × R50	(8095) V04 0216 029 237




Mtwo® Paper Points

STERILE

Suitable for the Mtwo® NiTi system

- Fully saturated points retain good material consistency
- Sterilised points in blister cells to support infection control
- 144 pieces: 4 papírové čepy v jednom blistru

	Size/Taper	29 mm
●	35/.04	(8096) V04 0225 429 035
●	40/.04	(8097) V04 0225 429 040
●	45/.04	(8098) V04 0225 429 045
●	50/.04	(8099) V04 0225 429 050
●	60/.04	(8100) V04 0225 429 060
●	30/.05	(8101) V04 0225 529 030
●	25/.06	(8102) V04 0225 629 025
●	30/.06	(8103) V04 0225 629 030
●	35/.06	(8104) V04 0225 629 035
●	40/.06	(8105) V04 0225 629 040
●	25/.07	(8106) V04 0225 729 025
	48 × 25/.06; 24 × 30/.05, 35/.04, 40/.04, 25/.07 each	(8107) V04 0225 029 457



pro•glider™

Nástroj je určený k vytvoření hladké cesty. Postačuje jediný nástroj velikosti 016.02 s variabilní konicitou.

- variabilní konicita zajišťuje trychtýřovitý vstup do kořenového kanálku a tím bezpečnější a pohodlnější přístup preparačních nástrojů
- nástroj je zhotovený z M-Wire slitiny (z níž jsou i nástroje Wave One a Protaper Next), což zaručuje vysokou flexibilitu a odolnost vůči cyklické únavě
- používá se při 300 ot./min. a točivém momentu 2 Ncm. Je možné ho kombinovat s libovolným reciprokačním nebo rotačním systémem na preparaci koř. kanálků
- je balený ve sterilních blistrech a koncipován pro jednorázové použití, což představuje následující výhody:
 - redukuje riziko zalomení nástroj na minimum
 - vyšší efektivita neopotřebovaného nástroje žádné riziko křížové kontaminaci
 - odpadá nutnost sterilizovat a organizovat nástroje



ProGlider® kónus .02 blister, baleno po 3 kusech

Size/Taper	21 mm	25 mm	31 mm
16/.02	(9130) 703 204 A0923 21 016	(9131) 703 204 A0923 25 016	(9132) 703 204 A0923 31 016




ProGlider® kónus .02 blister, baleno po 6 kusech

Size/Taper	21 mm	25 mm	31 mm
16/.02	(9133) 703 204 A0926 21 016	(9134) 703 204 A0926 25 016	(9135) 703 204 A0926 31 016



PathFile®

- rotační nástroje pro vytvoření hladké cesty "glide path"
- lépe zachovávají přirozenou anatomii kanálku
- více tolerují chyby v pracovním postupu jako je nesprávná délka (žádná apikální transpozice, pokud je pracovní doba příliš dlouhá a žádný schůdek, pokud je pracovní doba příliš krátká)
- redukují pracovní čas
- umožňují lékařům, kteří se nesespecializují na endodoncii, dosáhnout stejných výsledků jako specialisté



Pathfile® kónus .02 blister, baleno po 6 kusech

Size/Taper	21 mm	25 mm	31 mm
13/.02	(9136) 703 204 A0015 21 013	(9139) 703 204 A0015 25 013	(9142) 703 204 A0015 31 013
16/.02	(9137) 703 204 A0015 21 016	(9140) 703 204 A0015 25 016	(9143) 703 204 A0015 31 016
19/.02	(9138) 703 204 A0015 21 019	(9141) 703 204 A0015 25 019	(9144) 703 204 A0015 31 019
sada 013 - 019 (2x)	(9145) 703 204 A0015 21 S1	(9146) 703 204 A0015 25 S1	(9147) 703 204 A0015 31 S1



wave•one® Gold Glider

Reciprokační kořenový nástroj pro Glide Path

WaveOne Gold Glider® kónus .02 blister, baleno po 3 kusech

Size/Taper	21 mm	25 mm	31 mm
15/.02	(9175) 703 204 A0503 21 015	(9176) 703 204 A0503 25 015	(9177) 703 204 A0503 31 015



WaveOne Gold Glider® kónus .02 blister, baleno po 6 kusech

Size/Taper	21 mm	25 mm	31 mm
15/.02	(9178) 703 204 A0506 21 015	(9179) 703 204 A0506 25 015	(9180) 703 204 A0503 31 015



wave•one® Gold

Strojové Ni-Ti endodontické nástroje firmy Maillefer, které zásadně zjednodušují proces endodontické preparace.



WaveOne® kónus .07, .06, .05 blister, baleno po 3 kusech

Size/Taper	21 mm	25 mm	31 mm
20/.07 small	(9148) 703 204 A0700S 21 020	(9149) 703 204 A0700S 25 020	(9150) 703 204 A0700S 31 020
25/.07 primary	(9151) 703 204 A0700P 21 025	(9152) 703 204 A0700P 25 025	(9153) 703 204 A0700P 31 020
35/.06 medium	(9154) 703 204 A0700M 21 035	(9155) 703 204 A0700M 25 035	(9156) 703 204 A0700M 31 035
45/.05 large	(9157) 703 204 A0700L 21 045	(9158) 703 204 A0700L 25 045	(9159) 703 204 A0700L 31 045


STERILE



WaveOne® kónus .07, .06, .05 blister, baleno po 6 kusech

Size/Taper	21 mm	25 mm	31 mm
20/.07 small	(9160) 703 204 A0706S 21 020	(9161) 703 204 A0706S 25 020	(9162) 703 204 A0706S 31 020
25/.07 primary	(9163) 703 204 A0706P 21 025	(9164) 703 204 A0706P 25 025	(9165) 703 204 A0706P 31 020
35/.06 medium	(9166) 703 204 A0706M 21 035	(9167) 703 204 A0706M 25 035	(9168) 703 204 A0706M 31 035
45/.05 large	(9169) 703 204 A0706L 21 045	(9170) 703 204 A0706L 25 045	(9171) 703 204 A0706L 31 045

STERILE



WaveOne® kónus .07, .06, .05 blister, baleno po 4 sada

Size/Taper	21 mm	25 mm	31 mm
sada, každý 1x	(9172) 703 204 A0700 21 S1	(9173) 703 204 A07000 25 S1	(9174) 703 204 A0700 31 S1



MICRO-OPENER

Micro-opener baleno po 3 kusech




	Size/Taper		
□	0,10mm/.04	č.1	(9317) 703 A0350 04 010
□	0,15 mm/.04	č.2	(9318) 703 A0350 04 015
□	0,10 mm/.06	č.3	(9319) 703 A0350 06 010



MICRO-DEBRIDER

Micro-debrider (ruční H-file) baleno po 3 kusech



	Size/Taper		
□	0,20 mm/.02	č.1	(9196) 703 A0353 02 020
□	0,30 mm/.02	č.2	(9318) 703 A0353 02 030
□	10 mm/.06	č.3	(9319) 703 A0350 10 006 03



PROTAPER®NEXT

Nová generace Protaperů určená pro endo-odborníky hledající všestranný, flexibilní systém pro ošetření většiny kořenových kanálků.

Unikátní pohyb a větší flexibilita pro ještě náročnější, zakřivenější a úzké kanálky.

Riziko poškození nástroje je výrazně nižší a ještě lépe zachovávají původní anatomii kanálku.

Kratší klinická sekvence = méně nástrojů = kratší pracovní čas.

Speciální čtvercový průřez kanálku umožňuje vytlačování detritu v průběhu opracování.

Materiál M-WiRE, ze kterého je nástroj vyroben zlepšuje jeho flexibilitu, snižuje únavu materiálu a zachovává účinnost řezání.

Opět nechybí papírové a gutaperčové čepy, které přesně odpovídají tvaru kořenových nástrojů.



Protaper NEXT blister, baleno po 3 kusech

Size/Taper	21 mm	25 mm	31 mm
17/.04 X1	(6706) 703 X1 A0803 21 020	(6707) 703 X1 A0803 25 020	(6708) 703 X1 A0803 31 020
25/.06 X2	(6709) 703 X2 A0803 21 025	(6710) 703 X2 A0803 25 025	(6711) 703 X2 A0803 31 025
30/.07 X3	(6712) 703 X3 A0803 21 030	(6713) 703 X3 A0803 25 030	(6714) 703 X3 A0803 31 030
40/.06 X4	(6715) 703 X4 A0803 21 040	(6716) 703 X4 A0803 25 040	(6717) 703 X4 A0803 31 040
50/.06 X5	(6718) 703 X5 A0803 21 050	(6719) 703 X5 A0803 25 050	(6720) 703 X5 A0803 31 050
sada X1-X3, každý 1x	(6743) 703 A0803 X1-X3 21 SET	(6744) 703 A0803 X1-X3 25 SET	(6745) 703 A0803 X1-X3 31 SET

STERILE



Protaper NEXT blister, baleno po 6 kusech

Size/Taper	21 mm	25 mm	31 mm
17/.04 X1	(6721) 703 X1 A0800 21 020	(6722) 703 X1 A0800 25 020	(6723) 703 X1 A0800 31 020
25/.06 X2	(6724) 703 X2 A0800 21 025	(6725) 703 X2 A0800 25 025	(6726) 703 X2 A0800 31 025
30/.07 X3	(6727) 703 X3 A0800 21 030	(6728) 703 X3 A0800 25 030	(6729) 703 X3 A0800 31 030
40/.06 X4	(6730) 703 X4 A0800 21 040	(6731) 703 X4 A0800 25 040	(6732) 703 X4 A0800 31 040
50/.06 X5	(6733) 703 X5 A0800 21 050	(6734) 703 X5 A0800 25 050	(6735) 703 X5 A0800 31 050
sada X1-X3, každý 2x	(6736) 703 X1-X3 21 SET	(6737) 703 X1-X3 25 SET	(6738) 703 X1-X3 31 SET

STERILE



PROTAPER® UNIVERSAL RETREATMENT – ROTARY

Speciální systém 3 nástrojů na reendodocii

Rotační NiTi nástroje na reendodocii zaručující efektivitu a bezpečnost

Speciálně vyvinuté pro reendodocii k odstranění výplňového materiálu před opětovným opracováním kořenového kanálku

Protaper universal blister, baleno po 6 kusech

Size/Taper	
30/.09	(6740) 703 D1 16 009 030
25/.08	(6739) 703 D2 18 008 025
20/.07	(6741) 703 D3 22 007 020
sada d1 - D3, každý 2x	(6742) 703 D1 D3 SET



PROTAPER® GOLD

Nové nástroje ProTaper Gold přímo navazují na léty osvědčený, všeobecně uznávaný systém ProTaper Universal. K výhodám variabilní konicity přidává úžasnou flexibilitu a individuální přizpůsobení nástroje tvaru kanálku. Nový použitý materiál - speciální tepelně upravená slitina bez tvarové paměti umožňuje volné ohnutí nástroje a zachování tvaru - což zaručuje dokonalou adaptaci nástroje v kanálku, hlavně u dokončovacích nástrojů v apikální části.

ProTaper Gold versus ProTaper Universal

- 2,6x vyšší odolnost vůči cyklické únavě
- Vyšší flexibilita v průměru o 24% (nástroje umožňuje dodržet přirozenou anatomii kořenového kanálku)
- Zkrácený dřív ze 13 mm na 11 mm (vylepšuje intraorální viditelnost)
- Vyšší flexibilita nástroje díky patentované metalurgii, GOLD technologie
- Stejná sekvence nástrojů, technika opracování a řezná schopnost, jako u ProTaper®



Protaper gold SX® kónus .04 blister, baleno po 6 kusech

Size/Taper	19 mm
19/.04	(6699) 703 SX A0410 19


STERILE



Protaper gold® blister, baleno po 6 kusech

Size/Taper	21 mm	25 mm	31 mm
18/.02 S1	(6700) 703 S1 A0410 21 010	(6701) 703 S1 A0410 25 010	(6702) 703 S1 A0410 31 010
20/.04 S2	(6703) 703 S2 A0410 21 010	(6704) 703 S2 A0410 25 010	(6705) 703 S2 A0410 31 010
20/.07 F1	(9181) 703 F1 A0410 21 020	(9182) 703 F1 A0410 25 020	(9183) 703 F1 A0410 31 020
25/.08 F2	(9184) 703 F2 A0410 21 025	(9185) 703 F2 A0410 25 025	(9186) 703 F2 A0410 31 025
30/.09 F3	(9187) 703 F3 A0410 21 030	(9188) 703 F3 A0410 25 030	(9189) 703 F3 A0410 31 030
40/.06 F4	(9190) 703 F4 A0410 21 040	(9191) 703 F4 A0410 25 040	(9192) 703 F4 A0410 31 040
50/.05 F5	(9193) 703 F5 A0410 21 050	(9194) 703 F5 A0410 25 050	(9195) 703 F5 A0410 31 050
sada SX-F3	(9326) 703 A410 21 SET	(9327) 703 A410 25 SET	(9328) 703 A410 31 SET


STERILE



wave•one® GuttaCore, obturátory

- WaveOne Gold GuttaCore, blister
- balení 5 kusů obturátory + 1 kus verifikátor / 25 kusů obturátory + 5 kusů verifikátory

Size/Taper	balení 6 kusů	balení 30 kusů
 20/.07 small	(9338) 703 042 WOG 006 020	(9339) 703 042 WOG 030 020
 25/.07 primary	(9340) 703 042 WOG 006 025	(9341) 703 042 WOG 030 025
 35/.06 medium	(9342) 703 042 WOG 006 035	(9343) 703 042 WOG 030 035
 45/.05 large	(9344) 703 042 WOG 006 045	(9345) 703 042 WOG 030 045



wave•one® gutaperčové čepy

- WaveOne Gold Conform Fit
- balení 60 kusů

Size/Taper	28 mm
 20/.07 small	(9320) 703 042 000 07 020
 25/.07 primary	(9321) 703 042 000 07 025
 35/.06 medium	(9322) 703 042 000 06 035
 45/.05 large	(9323) 703 042 000 05 045
sada	(9324) 703 042 000 AS1




PROTAPER•NEXT® GUTTARORE, obturátory

- GuttaCore, obturátory pro systém Propeper Next
- balení 6 / 30 kusů







GuttaCore blister, baleno po 6 / 30 kusech

Size/Taper	balení 6 kusů	balení 30 kusů
 25/.06 X2	(9331) 703 042 X2 006 025	(9332) 703 042 X2 030 025
 30/.07 X3	(9333) 703 042 X3 006 030	(9334) 703 042 X3 030 030
 40/.06 X4	(9335) 703 042 X4 006 040	(9336) 703 042 X4 030 040
 50/.06 X5	(9337) 703 042 X5 006 050	



PROTAPER•NEXT® gutaperčové čepy

- ProTaper Next Conform Fit
- balení 60 kusů

Size/Taper	28 mm
 25/.06 X2	(9346) 703 042 X2 060 025
 30/.07 X3	(9347) 703 042 X3 060 030
 25/.06 /  30/.07 X2/X3	(9348) 703 042 X2/X3 060 S1
 50/.06 /  40/.06 X4/X5	(9349) 703 042 X4/X5 060 S1





Rotadent
ROTATING INSTRUMENTS
www.rotadent.eu


info@rotadent.eu


+420 388 323 991

Rotadent spol. s r.o.

**Špidrova 104
385 01 Vimperk**